



### CCRI FACS Core Unit instrument configurations

laser line	filter		
	Aria	Fortessa	Symphony
UV 355 nm		379/28	379/28
			450/50
		515/30	515/30
			580/20
			605/20
			670/25
		740/35	735/30
Violet 405 nm			810/40
	450/40	450/40	431/28
	525/50	525/50	525/50
			586/15
	610/20	605/12	605/40
	660/20	660/20	677/20
	710/50	710/50	710/50
Blue 488 nm	780/60	780/60	750/30
			810/40
	530/30	530/30	530/30
			610/20
	695/40	670/30	670/30
	710/50	710/50	710/50
			810/40
Yellow-green 561 nm	582/15	582/15	586/15
	610/20	610/20	610/20
	710/50		670/30
	780/60	780/60	780/60
Red 640 nm	670/30	670/14	670/30
	730/45	730/45	730/45
	780/60	780/60	780/60